



DEXXON<sup>®</sup>  
MEDICAL

DISPOSABLE POWDER FREE  
**GLOVE**

PERSONAL PROTECTIVE EQUIPMENT  
**A PROTECTION AND SAFETY SUPPLIER**

[www.dexxonmedical.com](http://www.dexxonmedical.com)





**RESTAURANTS**



**HOSPITALS**



**HEALTHCARE INSTITUTIONS**



**SHOPPING MARKETS**



**HOME TASKS**



**URBAN LIFE**



**OFFICES**



**CABINE CREWS**



**SCHOOLS**

The DEXXON MEDICAL GLOVES Termaplus Poly gloves prevent any allergic reactions among workers and customers. Made from polyethylene (LDPE) with 7.5/8 Micron, the gloves do not release any unwanted chemicals and therefore allow direct contact with food or any other sensitive items. As an excellent alternative to expensive gloves made from various other materials, DEXXON MEDICAL GLOVES are suitable for safety use in Food & Beverage industries, Consumer Services as SPA, beauty centers, barbers and daily life.



- ✓ Excellent fit, maximum sensitivity and dexterity for fine tasks. The entire glove is finely textured (including palms and fingers) to provide consistent grip on glassware, small objects, instruments and tools.
- ✓ Isolate odors and dust to keep your hands clean
- ✓ Ideal for all kinds of uses that require a safe and hygienic environment
- ✓ Ideal for almost any application including tattooing, food preparation, painting, cleaning, pet care, home improvement, hobbies, arts and crafts
- ✓ Convenient dispense pack keeps gloves clean and organized while providing easy access



## TECHNICAL DATA SHEET

### Artical Number

5001(S), 5002 (M ), 5003 (L), 5004 (XL)

### Packaging

100 PCS/Box | 24 Boxes/Carton | 2400 PCS/Carton

### Glove Specifications

Material: TPE | Color: Blue

Surface: Texturierte Fingerstip | Cuff Design: Rolled Edge | Latex: Free

### Thickness

Fingertips: 0,7 mm

### Physical Properties

EN 455-1 Freedom From Holes: 1,5 AQL

EN 455-2 Physical Properties : Long, 240 mm | Width, 110 mm

Dimensions : Force and break: 6,5 N

Force and break after : 6,2 N

EN 455-3 Powder Fee : 0,3 mg Per Glove

### Chemical Resistance for Degradation

TS EN 374-4 + EN 374-1

### Resistance to Permatation by Chemicals

EN 374-4, TS EN 16523-1, TS EN 659

### PPE Certification

Module A

### Virus Penetration

ISO16604

### Color and Sizes

**BLUE** SMALL | MEDIUM | LARGE | X-LARGE

### Dimensions

Box: B 196 mm x 110 mm x 48 mm Carton: B 420 mm x 235 mm x 310 mm

### Directive

Prepared according to 2016/425 Personal Protective Equipment Directive

### Expiration Time

36 Months after production

### Country of Manufacture

TURKEY

AVAILABLE COLOURS



AVAILABLE SIZES



100 GLOVES



24 BOXES




72 CARTONS



33 PALLETS

**1Box 100 PCS | 1Carton 24 Box | 1Palatte 72 Carton | 1Truck 57.024 Box**







DISPOSABLE POWDER FREE  
**GLOVE**

If you have any request,  
please don't hesitate to contact with me.

İstanbul Vizyon Park Ofis Blokları  
Yenibosna Merkez Mah. 29 Ekim Cad.  
Plaza:1 Kat:8 No:84  
Bahçelievler / İstanbul / Türkiye

[info@dexxon.com](mailto:info@dexxon.com)

[www.dexxon.com](http://www.dexxon.com)

[www.dexxonmedical.com](http://www.dexxonmedical.com)

[www.dexxonmedical.com](http://www.dexxonmedical.com)