



Disposable Mob Cap

FDA CE MADE
IN TURKEY

Product Details

Documents: CE, ISO 13485

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSSMC-01

Certificate Nr: NVA-EC-21022301

Product Name: Disposable surgical mob cap

Use: Disposable surgical non-woven Bouffant Cap

Material: % 100 Spunbond PP non-woven bonnet

Weight: 10-20 g/m²

Size: Standard size

Colour: Medical blue, white ,green

Classification: Class I



Disposable Mob Cap

FDA CE MADE
IN TURKEY

Performance indicators: Technical requirements for disposable mob cap
EN 13034:2005 +A1:2009,

Feature:

Cap with an elastic rubber band

Made of airy % 100 spunbond or SMS polypropylene non-woven

Excellent breathability and wearing comfort

Single-use, non-sterile

Antistatic

They are designed to provide additional strength while still offering breathability and comfort

Synthetic elastic offers a personalized fit.

Type : Surgical Supplies

Style : Single or double elastic

Application:

Hospital, medical, surgical, clean room, laboratory, beauty & food industry, catering industry

Etc

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).