



## Disposable Mob Cap

FDA CE MADE  
IN TURKEY

### Product Details

**Documents:** CE, ISO 13485

**Brand Name:** Dexxon Medical

**Model Nr:** DXNMD-DNSSMC-01

**Certificate Nr:** NVA-EC-21022301

**Product Name:** Disposable surgical mob cap

**Use:** Disposable surgical non-woven Bouffant Cap

**Material:** % 100 Spunbond PP non-woven bonnet

**Weight:** 10-20 g/m<sup>2</sup>

**Size:** Standard size

**Colour:** Medical blue, white ,green

**Classification:** Class I



## Disposable Mob Cap

FDA CE MADE  
IN TURKEY

**Performance indicators:** Technical requirements for disposable mob cap  
EN 13034:2005 +A1:2009,

**Feature:**

Cap with an elastic rubber band

Made of airy % 100 spunbond or SMS polypropylene non-woven

Excellent breathability and wearing comfort

Single-use, non-sterile

Antistatic

They are designed to provide additional strength while still offering breathability and comfort

Synthetic elastic offers a personalized fit.

**Type :** Surgical Supplies

**Style :** Single or double elastic

**Application:**

Hospital, medical, surgical, clean room, laboratory, beauty & food industry, catering industry

Etc

**Precautions:**

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).