



Disposable Coverall Type 4-B/ 5-B/ 6-B



MADE
IN TURKEY

Product Details

Product Description : Disposable Coverall PP

Brand Name: Dexxon Medical

Model Nr: DXNMD- DNSC-11

Certificated : CE, FDA and ISO13485

Product Name: Disposable Coverall Type 4B - 5B - 6B

Material quality: PE (polyethylene) film-coated non-woven fabric

Fabric Type: SMS Fabric microporous polypropylene + polyethylene

Usage: Standard protective coverall

Sterilization: Non-sterile

Weight: 65 - 67 g / m² (*It is the sum of the weight of the fabrics produced in m²)

Classification: Cat III, Type 4/5/6

Certifications :

Tested in accordance with PPE standard Type 4/5/6 Microscopus protective coverall with hood SMS polypropylene base material with PE microporous film coating

Related Directive : Personal Protective Equipment 2016/25 (EU)

Harmonised Standarts:

TYPE 4: EN 14605:2005+A1:2009

TYPE 5: EN 13982-1:2004+A1:2010

TYPE 6: EN 13034:2005+A1:2009

EN 13688:2013

EN 13034:2005 /A1:2010 (Type 6)

EN 14126 :2003 (Type 5-B, Type 6-B)

EN 1149-5:2008

EN 1073-2:2002

Regular Color: White

Basic Info

Name: Protective Coveralls / Disposable Coverall

Color : White

Size : S – M – L – XL – XXL – XXXL

Regular Material: PP/SMS



Disposable Coverall Type 4-B/ 5-B/ 6-B



MADE
IN TURKEY

Product features: Non-woven Non-Sterile Disposable Medical Protective Coveralls The product's resistance to synthetic blood penetration, breaking strength, and moisture permeability are above the national standard and have excellent performance for both safety and comfort.

Disposable Coverall PP

Disposable and Breathable laminated fabric with dust repellence at 1 micron level

Antistatic

Chemical resistant

Microporous polyethylene laminate

Comfort and durability guaranteed

Main features: medical-grade, bacteria-blocking, strong barrier; one-way protection and infection prevention.

Style/Pattern:

Taped zip for extra security

Easy fit hood for face masks

Elastic cuff

Product Specifications

Protective coveralls consists of a body piece, long sleeves, long trouser legs and a hood, attached a self-adhesive zipper cover strip.

Breathable material to help reduce heat build-up and promote comfortable wear. Elasticated waist and ankle for convenience

and freedom of movement for extra comfort.

Application / Typical areas of use

Protective coveralls designed to protect workers from hazardous substances or sensitive products and processes from contamination by people.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).