



Disposable Surgical Gown SMMS, Sterile

FDA CE MADE
IN TURKEY

Product Details

Product Description : SMMS Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSMG-SMS07

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21022305

Product Name: Disposable non-sterile surgical gown

Type : SMS Surgical Gown

Usage: Standard surgical gown

Sterilization: Non-sterile

Weight: 35-45 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS



Disposable Surgical Gown SMMS, Sterile

FDA CE MADE
IN TURKEY

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level : Level 1

Feature:

SMS Fabric
Liquid penetration
Non-Sterile

Application:

Our surgical gown is to be worn by medical staff during surgery.
The fabric feature is blue SMS fabric.
The apron features a loop neck closure for a custom fit without gaps.
Cuffs at the wrists for added protection.
Long sleeved gown with open back, under knee length
used in surgical/operating room



**Disposable Surgical
Gown SMS PE Laminated Reinforced,
Non Sterile**

  **MADE
IN TURKEY**

Product Details

Product Description : SMS PE Laminated Reinforced Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD- DNSLRSG-09

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21022307

Product Name: Disposable SMS Laminated Reinforced Sterile Surgical Gown

Material quality: PE (polyethylene) film-coated non-woven fabric

Type : Laminated Reinforced Surgical Gown Level 3 AAMI

Usage: Reinforced Surgical Gown

Sterilization: Non-sterile

Weight: 45 -55 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMS PE Laminated Reinforced Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS