



Disposable Surgical Gown SMS, Sterile

FDA CE MADE
IN TURKEY

Product Details

Product Description : SMS Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSMG-SSG18

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21061108

Product Name: Disposable non-sterile surgical gown

Type : SMS Surgical Gown

Usage: Standard surgical gown

Sterilization: Sterile by ethylene oxide

Weight: 35-45 g / m²

Classification: Class IS

Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMS



Disposable Surgical Gown SMS, Sterile

FDA CE MADE
IN TURKEY

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012
AAMI liquid barrier level : Level 1

Feature:

SMS Fabric
Liquid penetration
EO-Sterilized

Application:

Our surgical gown is to be worn by medical staff during surgery.
The fabric feature is blue SMS fabric.
The apron features a loop neck closure for a custom fit without gaps.
Cuffs at the wrists for added protection.
Long sleeved gown with open back, under knee length
used in surgical/operating room

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).