



Disposable Surgical Gown SS, Non-Sterile

  **MADE
IN TURKEY**

Product Details

Product Description : SS Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSMG-SS06

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21022304

Product Name: Disposable non-sterile surgical gown

Type : SS Surgical Gown

Usage: Standard surgical gown

Sterilization: Non-sterile

Weight: 35-45 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SS Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SS



Disposable Surgical Gown SS, Non-Sterile

  MADE
IN TURKEY

Performance indicators:

Technical requirements for disposable surgical gown
EN 13795 : 2019

Feature:

Disposable
SS Fabric
Non-Sterile

Application:

Our surgical gown is designed to be worn by healthcare professionals during surgical procedures.

Fabric feature is blue SS fabric.

The apron features a loop neck closure for an easily adjustable, gapless custom fit.

Cuffs at the wrists for added protection

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).