



DEXXON
MEDICAL

PRODUCT CATALOG

High quality standards and professionalism in
medical textile products





DEXXON MEDICAL

ABOUT US

We are proud to be the choice of global companies thanks to our qualified service and business partnership that we provide to our customers with our years of experience and knowledge.

DEXXON MEDICAL brand , is one of leading Medical Garment exporters ,combining high technology and experience in its products.

Our company incorporates a rich medical textile culture with its global customer portfolio and exports to Europe, Asia and North America.

WHO ARE WE ?

With years of experience and knowledge in the finished textile industry, we pride ourselves on being the choice of global companies thanks to the high quality service and business understanding we offer our customers in the medical textile industry.

DEXXON MEDICAL is active in the health sector and has combined its experience and know-how in the field of ready-made goods with its innovative management approach and its strong supplier network.

The company is one of the world's leading companies in the field of sterile single-use medical devices with 100% customer satisfaction. All production facilities of our company in accordance with the rules of the principles of the quality management system ISO 13485 established in the medical device industry; It is audited and certified by international organizations with regard to social compliance and codes of conduct, health and safety, and ethical codes.

FFP2 NR Face Masks

Product details

- Brand Name : Dexxon Medical
- Model Nr : DXNMD-NRFM05
- Product Name : Disposable face mask
- Use : Disposable face mask
- Material : Spunbond/Cotton/Melt-Blown
- Size : Spunbond/Melt-Blown/Spunbond
- Color : Standard
- Packaging : White (all colors possible)
- Classification : 50 pcs/package

*All test and inspection results are available and Can be sent upon request.

FFP2 (10 Adet)



FFP2 (25 Adet)



FFP2 (50 Adet)



TURUNCU FFP2 MASKE



PEMBE FFP2 MASKE



KAHVERENGİ FFP2 MASKE



SIYAH FFP2 MASKE



BEYAZ FFP2 MASKE



GRİ FFP2 MASKE



KIRMIZI FFP2 MASKE



BORDO FFP2 MASKE



YEŞİL FFP2 MASKE



FFP3 NR Face Masks

Product details

- Documents : CE, ISO
- Brand Name : Dexxon Medical
- Model Nr : DXNMD-NRFM04
- Product Name : Disposable face mask
- Use : Disposable face mask
- Material : Spunbond/Cotton/Melt-Blown
Melt-Blown/Spunbond
- Size : Standard
- Color : White (all colors possible)
- Packaging : 50 pcs/package
- Classification : FFP2

*All test and inspection results are available and Can be sent upon request.

FFP3 (10 Adet)



FFP3 (25 Adet)



FFP3 (50 Adet)





Disposable Surgical Mask Type IIR



Type IIR - EN14683, the technical data are:

- Splash-proof layer protection against blood and body fluids
- High filter capacity
- Easy breathing
- The adjustable nose wire, flexible and soft, fits the face perfectly
- Mask with elastic rubber
- Breathable, makes breathing easier
- Anti-pilling, 3 layers and 3 layers / wrinkles
- Does not contain latex
- Sterile or non-sterile
- Disposable
- Each product is individually packaged
- Max box content 50 pcs

A-) OUTER LAYER: Non-woven antibacterial fabric.
 Splash resistant to EN 14683:2019 standards.
 30 g/m² Polypropylene Spunbond non -woven fabric

B-) MIDDLE LAYER: High barrier melt blown filter
 BFE>99 %

25 g/m² Polypropylene Meltblown non-woven fabric

C-) INNER LAYER: Sensitive to skin and moisture absorbing quality

20 g/m² Polypropylene Spunbond non -woven fabric

Medical face mask - Type IIR - Medium splash protection

Protecting health care workers and patients from infectious diseases



Disposable Surgical Mask Type IIR



Main Features and Benefits

- More comfort and functionality,
- Provides softness and sweat absorption
- Provides extra space between the lips and the mask for better breathing and is skinfriendly
- Reduces leakage and offers a higher filter capacity
- Splash-proof layer
- Developed for surgery and use in situations where there is a risk of splashing
- Also ideal for dental applications

Recommended for the following

- Improving infection control
- Prevention of cross contamination risk
- Medical face masks prevent large particles (e.g. saliva, mucus) consumed by the user from reaching the patient or the work environment

DISPOSABLE PROTECTIVE GOWN





Product Details

Product Description : SS Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSMG-SS06
Certificated : CE and ISO13485
Certificate Nr: NVA-EC-21022304
Product Name: Disposable non-sterile surgical gown
Type : SS Surgical Gown
Usage: Standard surgical gown
Sterilization: Non-sterile
Weight: 35-45 g / m2
Classification: Class I
Regular Color: Medical Blue

Basic Info

Name: SS Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SS

Performance indicators:

Technical requirements for disposable surgical gown
 EN 13795 : 2019

Feature:

Disposable
 SS Fabric
 Non-Sterile

Application:

Our surgical gown is designed to be worn by healthcare professionals during surgical procedures.
 Fabric feature is blue SS fabric.
 The apron features a loop neck closure for an easily adjustable, gapless custom fit.
 Cuffs at the wrists for added protection

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SS Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSMG-SS04
Certificated : CE and ISO13485
Certificate Nr: NVA-EC-21052601
Product Name: Disposable sterile surgical gown
Type : SS Surgical Gown
Usage: Standard surgical gown
Sterilization: Sterile by ethylene oxide
Weight: 35-45 g / m2
Classification: Class IS
Regular Color: Medical Blue

Basic Info

Name: SS Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SS

Performance indicators:

Technical requirements for disposable surgical gown
 EN 13795 : 2019

Feature:

Disposable
 SS Fabric
 EO-Sterilized

Application:

Our surgical gown is designed to be worn by healthcare professionals during surgical procedures.
 Fabric feature is blue SS fabric.
 The apron features a loop neck closure for an easily adjustable, gapless custom fit.
 Cuffs at the wrists for added protection

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SMS Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSMG-SMS07
Certificated : CE and ISO13485
Certificate Nr: NVA-EC-21022305
Product Name: Disposable non-sterile surgical gown
Type : SMS Surgical Gown
Usage: Standard surgical gown
Sterilization: Non-sterile
Weight: 35-45 g / m2
Classification: Class I
Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SMMS

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012
AAMI liquid barrier level : Level 1

Feature:

SMS Fabric
 Liquid penetration
 Non-Sterile

Application:

Our surgical gown is to be worn by medical staff during surgery. The fabric feature is blue SMS fabric. The apron features a loop neck closure for a custom fit without gaps. Cuffs at the wrists for added protection. Long sleeved gown with open back, under knee length used in surgical/operating room



Product Details

Product Description : SMS Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSMG-SSG18
Certificated : CE and ISO13485
Certificate Nr: NVA-EC-21061108
Product Name: Disposable non-sterile surgical gown
Type : SMS Surgical Gown
Usage: Standard surgical gown
Sterilization: Sterile by ethylene oxide
Weight: 35-45 g / m2
Classification: Class IS
Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SMS

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012
AAMI liquid barrier level : Level 1

Feature:
 SMS Fabric
 Liquid penetration
 EO-Sterilized

Application:
 Our surgical gown is to be worn by medical staff during surgery.
 The fabric feature is blue SMS fabric.
 The apron features a loop neck closure for a custom fit without gaps.
 Cuffs at the wrists for added protection.
 Long sleeved gown with open back, under knee length
 used in surgical/operating room

Precautions:
 The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SMMS Surgical Gown / Scrub Suits

Brand Name: Dexxon Medical

Model Nr: DXNMD-DN55S14

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21061104

Product Name: Disposable non-sterile short-sleeved patient gown

Type : SMMS Surgical Gown / Scrub Suits

Usage: Standard patient gown

Sterilization: Non-sterile

Weight: 35-40 g / m2

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMMS Surgical Gown / Scrub Suits

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level : Level 1

Feature:

Disposable and Breathable

Spunbond (SMMS) consists of 100% PP non-woven fabric

Economical disposable short sleeves

Antistatic

Each product is packed individually

EN 13795

Style/Pattern: Two-piece set including shirt and pants.

Shirt:

Collar: Around collar or V-collar

Sleeves: Without or short sleeves

Pocket: Available

Pants: With tie on the waist, pockets are available.

Application: Hospital, laboratory etc

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SMS PE Laminated Surgical Isolation Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSLIG13
Certificated : CE and ISO13485
Certificate Nr:
Product Name: Disposable Isolation Clothing (Coveralls)
Type : Type 6-B Offering limited protective performance against liquid chemicals
Usage: Surgical Isolation Gown
Sterilization: Sterile
Weight: 45 -55 g / m2
Classification: CAT III
Regular Color: Medical Blue

Basic Info

Name: SMS PE Laminated Surgical Isolation Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SS or SMS

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products
 EN 14126:2003 AC:2004 , EN 13034:2005 +A1:2009, EN 13795-1:2019; AAMI PB70: 2012

Type : Type 6-B

Level: Level 3

Feature:

Disposable laminated surgical isolation gown

Scope of Application:It is used for general isolation in clinics, wards, and laboratories of medical institutions. It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection

Structural Composition:The product is mainly made of non-woven cloth through cutting and sewing. It is provided in a non-sterile condition and only for single use.

Clothing Features and Overview:

The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure. Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure.

It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.

Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Comfortable fabric: New materials, soft, light, breathable.

Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.

Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.

Characteristics

Excellent water repellence to protect surgeons from contamination.

Prevention of cross contamination by reducing lint

Resistant to liquids with its barrier feature

Resistant to moisture and bacterial penetration

Puncture, tear and abrasion resistant.

Free from toxic substances

EO-Sterilized

Application:

This gown can be used during the surgery operations.

This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%)



**Disposable Surgical Isolation Gown
SMS PE Laminated
Type 6-B Non-Sterile**

FDA CE MADE IN TURKEY

Product Details

- Product Description :** SMS PE Laminated Surgical Isolation Gown
- Brand Name:** Dexxon Medical
- Model Nr:** DXNMD-DNSLIG13
- Certificated :** CE and ISO13485
- Certificate Nr:** NVA-EC- 21061103
- Product Name:** Disposable Isolation Clothing (Coveralls)
- Type :** Type 6-B Offering limited protective performance against liquid chemicals
- Usage:** Surgical Isolation Gown
- Sterilization:** Non-Sterile
- Weight:** 45 -55 g / m2
- Classification:** CAT III
- Regular Color:** Medical Blue

Basic Info

- Name:** SMS PE Laminated Surgical Isolation Gown
- Color :** Blue
- Size :** S.M.L.XL.XXL.XXXL
- Regular Material:** PP/SS or SMS



**Disposable Surgical Isolation Gown
SMS PE Laminated
Type 6-B Non-Sterile**

FDA CE MADE IN TURKEY

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products
EN 14126:2003 AC:2004 , EN 13034:2005 +A1:2009, EN 13795-1:2019; AAMI PB70: 2012

Type : Type 6-B

Level: Level 3

Feature:

Disposable laminated surgical isolation gown

Scope of Application:It is used for general isolation in clinics, wards, and laboratories of medical institutions. It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection

Structural Composition:The product is mainly made of non-woven cloth through cutting and sewing. It is provided in a non-sterile condition and only for single use.

Clothing Features and Overview:

The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure. Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure.

It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.

Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Comfortable fabric: New materials, soft, light, breathable.

Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.

Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.

Characteristics

Excellent water repellence to protect surgeons from contamination.

Prevention of cross contamination by reducing lint

Resistant to liquids with its barrier feature

Resistant to moisture and bacterial penetration

Puncture, tear and abrasion resistant.

Free from toxic substances

Non-Sterile

Application:

This gown can be used during the surgery operations.

This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SMS PE Laminated Reinforced Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD- DNSLRSG-09
Certificated : CE and ISO13485
Certificate Nr: NVA-EC-21022307
Product Name: Disposable SMS Laminated Reinforced Sterile Surgical Gown
Material quality: PE (polyethylene) film-coated non-woven fabric
Type : Laminated Reinforced Surgical Gown Level 3 AAMI
Usage: Reinforced Surgical Gown
Sterilization: Non-sterile
Weight: 45 -55 g / m2
Classification: Class I
Regular Color: Medical Blue

Basic Info

Name: SMS PE Laminated Reinforced Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SMMS

Performance indicators: Technical requirements for disposable polyethylene laminated reinforced isolation clothing products / EN 13795-1:2019; AAMI PB70: 2012
Level : Level 3

Feature:

SMS Non-Woven Fabric
 Excellent water repellence to protect surgeons from contamination.
 Prevention of cross contamination by reducing lint
 Resistant to liquids with its barrier feature
 Resistant to moisture and bacterial penetration
 Puncture, tear and abrasion resistant.
 Free from toxic substances
 Non-Sterile

Application:

This gown can be used during the surgery operations.
 This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SMS PE Laminated Reinforced Surgical Gown
Brand Name: Dexion Medical
Model Nr: DXNMD- DNSLRSG-09
Certificated : CE and ISO13485
Certificate Nr:
Product Name: Disposable SMS Laminated Reinforced Sterile Surgical Gown
Material quality: PE (polyethylene) film-coated non-woven fabric
Type : Laminated Reinforced Surgical Gown Level 3 AAMI
Usage: Reinforced Surgical Gown
Sterilization: Sterile by ethylene oxide
Weight: 45 -55 g / m2
Classification: Class I
Regular Color: Medical Blue

Basic Info

Name: SMS PE Laminated Reinforced Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SMMS

Performance indicators: Technical requirements for disposable polyethylene laminated reinforced isolation clothing products / EN 13795-1:2019; AAMI PB70: 2012
Level : Level 3

Feature:

SMS Non-Woven Fabric
 Excellent water repellence to protect surgeons from contamination.
 Prevention of cross contamination by reducing lint
 Resistant to liquids with its barrier feature
 Resistant to moisture and bacterial penetration
 Puncture, tear and abrasion resistant.
 Free from toxic substances
 EO-Sterilized

Application:

This gown can be used during the surgery operations.
 This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SMMS Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSMG-SMS07
Certificated : CE and ISO13485
Certificate Nr: NVA-EC-21022305
Product Name: Disposable non-sterile surgical gown
Type : SMS Surgical Gown
Usage: Standard surgical gown
Sterilization: Non-sterile
Weight: 35-45 g / m2
Classification: Class I
Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SMMS

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012
AAMI liquid barrier level : Level 1

Feature:

SMS Fabric
 Liquid penetration
 Non-Sterile

Application:

Our surgical gown is to be worn by medical staff during surgery.
 The fabric feature is blue SMS fabric.
 The apron features a loop neck closure for a custom fit without gaps.
 Cuffs at the wrists for added protection.
 Long sleeved gown with open back, under knee length
 used in surgical/operating room



Disposable Surgical Gown SMMS, Non-Sterile

Product Details

Product Description : SMMS Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSMG-SMMS08

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21050501

Product Name: Disposable non-sterile surgical gown

Type : SMMS Surgical Gown

Usage: Standard surgical gown

Sterilization: Non-sterile

Weight: 35-45 g / m2

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMMS Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS



Disposable Surgical Gown SMMS, Non-Sterile

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level: Level 2

Feature:

Disposable surgical isolation gown-standard performance

SMMS Fabric

Liquid penetration

Scope of Application: It is used for general isolation in clinics, wards, and laboratories of medical institutions. It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection

Structural Composition: The product is mainly made of non-woven cloth through cutting and sewing. It is provided in a non-sterile condition and only for single use.

Clothing Features and Overview:

The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure.

Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure.

It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.

Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Characteristics

Comfortable fabric: New materials, soft, light, breathable.

Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.

Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.

Application:

Our surgical gown is to be worn by medical staff during surgery.

The fabric feature is blue SMMS fabric.

The exam gowns are durable, soft, light, comfortable, abrasion and static electricity resistant.

The exam gowns are air and vapor permeable, while resisting water fluid and aerosol.

The gown is made from 100% high density non woven, containing no additives, adhesives.

The apron features a loop neck closure for a custom fit without gaps.

Cuffs at the wrists for added protection.

All-in-one, front opening

Long sleeved gown with open back, under knee length

Used in surgical/operating room

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



DEXXON MEDICAL

**SMMS Surgical Gown,
Disposable Patient Use
SMMS Scrub Suit, Non-Sterile**

FDA CE MADE IN TURKEY

Product Details

Product Description : SMS Surgical Gown / Scrub Suits

Brand Name: Dexxon Medical

Model Nr: DXNMD-DN55S14

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21061104

Product Name: Disposable non-sterile short-sleeved patient gown

Type : SMS Surgical Gown / Scrub Suits

Usage: Standard patient gown

Sterilization: Non-sterile

Weight: 35-40 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown / Scrub Suits

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS



DEXXON MEDICAL

**SMMS Surgical Gown,
Disposable Patient Use
SMMS Scrub Suit, Non-Sterile**

FDA CE MADE IN TURKEY

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level : Level 1

Feature:

Disposable and Breathable

Spunbond (SMS) consists of 100% PP non-woven fabric

Economical disposable short sleeves

Antistatic

Each product is packed individually

EN 13795

Style/Pattern: Two-piece set including shirt and pants.

Shirt:

Collar: Around collar or V-collar

Sleeves: Without or short sleeves

Pocket: Available

Pants: With tie on the waist, pockets are available.

Application: Hospital, laboratory etc

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



**Disposable Surgical Gown
SMMS PE Laminated Reinforced,
Sterile**

Product Details

Product Description : SMMS PE Laminated Reinforced Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSMG- DSLRSG12
Certificated : CE and ISO13485
Certificate Nr: NVA-EC- 21052604
Product Name: Disposable SMMS Laminated Reinforced Sterile Surgical Gown
Type : Laminated Reinforced Surgical Gown Level 3 AAMI
Usage: Reinforced Surgical Gown
Sterilization: Sterile by ethylene oxide
Weight: 45 -55 g / m2
Classification: Class IS
Regular Color: Medical Blue

Basic Info

Name: SMMS PE Laminated Reinforced Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: PP/SMMS



**Disposable Surgical Gown
SMMS PE Laminated Reinforced,
Sterile**

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012
AAMI liquid barrier level : Level 3
Scope of Application: It is used for general isolation in clinics, wards, and laboratories of medical institutions. It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection
Structural Composition: The product is mainly made of non-woven cloth through cutting and sewing. It is provided in a non-sterile condition and only for single use.
Clothing Features and Overview:
 The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure. Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure. It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.
Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Characteristics
Comfortable fabric: New materials, soft, light, breathable.
Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.
Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.
Feature:
 Disposable laminated reinforced surgical gown performance
 SMMS Non-Woven Fabric
 Excellent water repellence to protect surgeons from contamination.
 Prevention of cross contamination by reducing lint
 Resistant to liquids with its barrier feature
 Resistant to moisture and bacterial penetration
 Puncture, tear and abrasion resistant.
 Free from toxic substances
 EO-Sterilized
Application:
 This gown can be used during the surgery operations.
 This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:
 The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Disposable Surgical Gown
SSMMS PE Laminated Extra Reinforced,
Non-Sterile

FDA CE MADE IN TURKEY

Product Details

Product Description : SSMMS PE Laminated Extra Reinforced Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSLERSG15

Certificated : CE and ISO13485

Certificate Nr: NVA-EC- 21061105

Product Name: Disposable SSMMS Laminated Reinforced Sterile Surgical Gown

Material quality: PE (polyethylene) film-coated non-woven fabric

Type : Laminated Reinforced Surgical Gown Level 3 AAMI ,Disposable laminated

Usage: Reinforced Surgical Gown

Sterilization: Non-Sterile

Weight: 55 g / m2

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SSMMS PE Laminated Reinforced Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: : PP/SSMMS



Disposable Surgical Gown
SSMMS PE Laminated Extra Reinforced,
Non-Sterile

FDA CE MADE IN TURKEY

SSMMS technology, used as the main material and as reinforcement in critical areas, provides high comfort and breathability for an effective gown.

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level : Level 3

Scope of Application:It is used for general isolation in clinics, wards, and laboratories of medical institutions. It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection

Structural Composition:The product is mainly made of non-woven cloth through cutting and sewing. It is provided in a non-sterile condition and only for single use.

Clothing Features and Overview:

The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure.

Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure.

It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.

Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Characteristics

Comfortable fabric: New materials, soft, light, breathable.

Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.

Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.

Excellent water repellence to protect surgeons from contamination.

Prevention of cross contamination by reducing lint

Resistant to liquids with its barrier feature

Resistant to moisture and bacterial penetration

Puncture, tear and abrasion resistant.

Free from toxic substances

Non-Sterile

Application:

This gown can be used during the surgery operations.

This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : SSMMS PE Laminated Extra Reinforced Surgical Gown
Brand Name: Dexxon Medical
Model Nr: DXNMD-DSLERSG16
Certificated : CE and ISO13485
Certificate Nr: NVA-EC- 21061106
Product Name: Disposable SSMMS Laminated Reinforced Sterile Surgical Gown
Material quality: PE (polyethylene) film-coated non-woven fabric
Type : Laminated Reinforced Surgical Gown Level 3 AAMI ,Disposable laminated extra reinforced surgical gown performance
Usage: Reinforced Surgical Gown
Sterilization: Sterile by ethylene oxide
Weight: 55 g / m2
Classification: Class IS
Regular Color: Medical Blue

Basic Info

Name: SSMMS PE Laminated Reinforced Surgical Gown
Color : Blue
Size : S.M.L.XL.XXL.XXXL
Regular Material: : PP/SSMMS

SSMMS technology, used as the main material and as reinforcement in critical areas, provides high comfort and breathability for an effective gown.

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level : Level 3

Scope of Application:It is used for general isolation in clinics, wards, and laboratories of medical institutions. It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection

Structural Composition:The product is mainly made of non-woven cloth through cutting and sewing. It is provided in a non-sterile condition and only for single use.

Clothing Features and Overview:
 The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure.
 Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure.
 It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.
Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Characteristics
Comfortable fabric: New materials, soft, light, breathable.
Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.
Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.
 Excellent water repellence to protect surgeons from contamination.
 Prevention of cross contamination by reducing lint
 Resistant to liquids with its barrier feature
 Resistant to moisture and bacterial penetration
 Puncture, tear and abrasion resistant.
 Free from toxic substances
 EO-Sterilized

Application:
 This gown can be used during the surgery operations.
 This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:
 The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).

DISPOSABLE COVERALL





Product Details

Product Description: Disposable Coverall PP
Brand Name: Dexxon Medical
Model Nr: DXNMD- DNSC-12
Certificated : CE, FDA and ISO13485
Product Name: Disposable Coverall Type 3B - 4B
Material quality: PE (polyethylene) film-coated non-woven fabric
Fabric Type: SMS Fabric microporous polypropylene + polyethylene
Usage: Standard protective coverall
Sterilization: Non-sterile
Weight: 55 - 57 g / m2 (*It is the sum of the weight of the fabrics produced in m2)
Classification: Cat III Type 3B-4B
 Type 3/4 Microscopu protective coverall with hood
 SMS polypropylene base material with PE microporous film coating
Certifications :
 Tested in accordance with PPE standard
Related Directive : Personal Protective Equipment 2016/25 (EU)
Harmonised Standarts:
 TYPE 3B/4B : EN 14605:2005+A1:2009
 EN 13688:2013
 EN 14126 :2003
 EN 1149-5:2008
 EN 1073-2:2002
Regular Color: White

Basic Info

Name: Protective Coveralls / Disposable Coverall
Color : White
Size : S – M – L – XL – XXL – XXXL
Regular Material: PP/SMS

Product features: Non-woven Non-Sterile Disposable Medical Protective Coveralls The product's resistance to synthetic blood penetration, breaking strength, and moisture permeability are above the national standard and have excellent performance for both safety and comfort.

Disposable Coverall PP
 Disposable and Breathable laminated fabric with dust repellence at 1 micron level
 Antistatic
 Chemical resistant
 Microporous polyethylene laminate
 Comfort and durability guaranteed
Main features: medical-grade, bacteria-blocking, strong barrier; one-way protection and infection prevention.

Style/Pattern:

Taped zip for extra security
 Easy fit hood for face masks
 Elastic cuff

Product Specifications

Protective coveralls consists of a body piece, long sleeves, long trouser legs and a hood, attached a self-adhesive zipper cover strip.
 Breathable material to help reduce heat build-up and promote comfortable wear.
 Elasticated waist and ankle for convenience and freedom of movement for extra comfort.

Application / Typical areas of use

Protective coveralls designed to protect workers from hazardous substances or sensitive products and processes from contamination by people.

Scope of application: For occupational safety protection of high-polluting production personnel, to prevent dust and other substances from splashing on people.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : Disposable Coverall PP

Brand Name: Dexxon Medical

Model Nr: DXNMD- DNSC-11

Certificated : CE, FDA and ISO13485

Product Name: Disposable Coverall Type 4B - 5B - 6B

Material quality: PE (polyethylene) film-coated non-woven fabric

Fabric Type: SMS Fabric microporous polypropylene + polyethylene

Usage: Standard protective coverall

Sterilization: Non-sterile

Weight: 65 - 67 g / m2 (*It is the sum of the weight of the fabrics produced in m2)

Classification: Cat III, Type 4/5/6

Certifications :

Tested in accordance with PPE standard Type 4/5/6 Microscop protective coverall with hood SMS polypropylene base material with PE microporous film coating

Related Directive : Personal Protective Equipment 2016/25 (EU)

Harmonised Standarts:

TYPE 4: EN 14605:2005+A1:2009

TYPE 5: EN 13982-1:2004+A1:2010

TYPE 6: EN 13034:2005+A1:2009

EN 13688:2013

EN 13034:2005 /A1:2010 (Type 6)

EN 14126 :2003 (Type 5-B, Type 6-B)

EN 1149-5:2008

EN 1073-2:2002

Regular Color: White

Basic Info

Name: Protective Coveralls / Disposable Coverall

Color : White

Size : S - M - L - XL - XXL - XXXL

Regular Material: PP/SMS

Product features: Non-woven Non-Sterile Disposable Medical Protective Coveralls The product's resistance to synthetic blood penetration, breaking strength, and moisture permeability are above the national standard and have excellent performance for both safety and comfort.

Disposable Coverall PP

Disposable and Breathable laminated fabric with dust repellence at 1 micron level

Antistatic

Chemical resistant

Microporous polyethylene laminate

Comfort and durability guaranteed

Main features: medical-grade, bacteria-blocking, strong barrier; one-way protection and infection prevention.

Style/Pattern:

Taped zip for extra security

Easy fit hood for face masks

Elastic cuff

Product Specifications

Protective coveralls consists of a body piece, long sleeves, long trouser legs and a hood, attached a self-adhesive zipper cover strip.

Breathable material to help reduce heat build-up and promote comfortable wear. Elasticated waist and ankle for convenience

and freedom of movement for extra comfort.

Application / Typical areas of use

Protective coveralls designed to protect workers from hazardous substances or sensitive products and processes from contamination by people.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of

direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



Product Details

Product Description : Disposable Coverall PP
Brand Name: Dexion Medical
Model Nr: DXNMD- DNSC-10
Certificated : CE, FDA and ISO13485
Product Name: Disposable Coverall Type 5B-6B
Material quality: PE (polyethylene) film-coated non-woven fabric
Fabric Type: SMS Fabric microporous polypropylene + polyethylene
Usage: Standard protective coverall
Sterilization: Non-sterile
Weight: 55 - 57 g / m2 (*It is the sum of the weight of the fabrics produced in m2)
Classification: Cat III
Certifications :
 Tested in accordance with PPE standard
 Cat III Type 5B/6B
 Type 5B/6B Microscopus protective coverall with hood
 SMS polypropylene base material with PE microporous film coating
 Related Directive : Personal Protective Equipment 2016/25 (EU)
Harmonised Standarts:
 TYPE 5: EN 13982-1:2004+A1:2010
 TYPE 6: EN 13034:2005+A1:2009
 EN 13688:2013
 EN 13034:2005 /A1:2010 (Type 6)
 EN 14126 :2003 (Type 5-B, Type 6-B)
 EN 1149-5:2008
Regular Color: White

Basic Info

Name: Protective Coveralls / Disposable Coverall
Color : White
Size : S – M – L – XL – XXL – XXXL
Regular Material: PP/SMS

Product features: Non-woven non-sterile disposable medical protective coveralls
 The product's resistance to synthetic blood penetration, breaking strength, and moisture permeability are above the national standard and have excellent performance for both safety and comfort.

- Disposable Coverall PP
- Disposable and Breathable
- Antistatic
- Chemical resistant
- Microporous polyethylene laminate
- Comfort and durability guaranteed

Style/Pattern:

- Concealed zip
- Double stitched with an additional overlap
- Elastic hem at ankles, wrists and waist
- The three-piece hood is specially designed to accommodate common respiratory protection materials
- Thumb or finger loop

Application: Hospital, laboratory etc

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).

Disposable Mob Cap





Product Details

Documents: CE, ISO 13485
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSSMC-01
Certificate Nr: NVA-EC-21022301
Product Name: Disposable surgical mob cap
Use: Disposable surgical non-woven Bouffant Cap
Material: % 100 Spunbond PP non-woven bonnet
Weight: 10-20 g/m2
Size: Standard size
Colour: Medical blue, white ,green
Classification: Class I

Performance indicators: Technical requirements for disposable mob cap EN 13034:2005 +A1:2009,

Feature:
 Cap with an elastic rubber band
 Made of airy % 100 spunbond or SMS polypropylene non-woven
 Excellent breathability and wearing comfort
 Single-use, non-sterile
 Antistatic
 They are designed to provide additional strength while still offering breathability and comfort
 Synthetic elastic offers a personalized fit.

Type : Surgical Supplies

Style : Single or double elastic

Application:

Hospital, medical, surgical, clean room, laboratory, beauty & food industry, catering industry Etc

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).

Disposable Shoe Cover

**Disposable PE
Laminated Shoe Cover**



**Disposable PE
Shoe Cover**



**Disposable
PP Shoe Cover**





Product Details

Documents: CE, ISO 13485
Brand Name: Dexxon Medical
Model Nr: DXNMD-DNSSC-02
Certificate Nr: NVA-EC-21022302
Product Name: Disposable PE Laminated shoe cover
Use: Disposable surgical shoe cover
Material: % 100 Spunbond PP non-woven shoe cover
Material quality: PE (polyethylene) film-coated non-woven fabric
Fabric Type: SMS fabric microporous polypropylene + polyethylene
Usage: Standard protective shoe cover
Sterilization: Non-sterile
Weight: 55 g / m2
Size: Standard size
Colour: Medical blue, white

Performance indicators: Technical requirements for disposable laminated shoe cover EN 14126 :2003 (Type 5-B, Type 6-B)
 Tested in accordance with PPE standard Cat III ; Type 5/6
 Type 5/6 Microscopus protective coverall with hood
 SMS polypropylene base material with PE microporous film coating
Related Directive : Personal Protective Equipment 2016/25 (EU)
Product features: Non-woven non-sterile disposable medical protective coveralls
 The product's resistance to synthetic blood penetration, breaking strength, and moisture permeability are above the national standard and have excellent performance for both safety and comfort.
 Disposable and Breathable
 Spunbond (SS or SMS) consists of 100% PP laminated non-woven fabric
 Antistatic
 Chemical resistant
 Microporous polyethylene laminate
 Comfort and durability guaranteed
Performance indicators: Technical requirements for disposable laminated shoe cover
Application: Hospital, medical, surgical, clean room,
Precautions:
 The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%



Disposable PE Shoe Cover

FDA CE MADE IN TURKEY

Product Details

Documents: CE, ISO 13485
Brand Name: Dexxon Medical
Model Nr: DXNMD-DPESC04
Certificate Nr:
Product Name: Disposable PE shoe cover
Use: Disposable shoe cover
Material: % 100 PE shoe cover
Material quality: Made of PE
Fabric Type: SS or SMS fabric
Sterilization: Non-sterile
Weight: 25 g/m²
Classification:
Size: Free Size
Colour: Blue



Disposable PE Shoe Cover

FDA CE MADE IN TURKEY

Product features:

Disposable
 100% PE fabric
 Antistatic
 Non-slip
 Dust-proof

Performance indicators:

Our products are mainly used in food factories, lab room, school, public places, medical treatment hygiene for preventing dust and static.



Disposable PP Shoe Cover



Disposable PP Shoe Cover

Product Details

Documents: CE, ISO 13485

Brand Name: Dexxon Medical

Model Nr: DXNMD-DPPSC-03

Certificate Nr:

Product Name: Disposable PP shoe cover

Use: Disposable shoe cover

Material: % 100 Spunbond PP non-woven shoe cover

Material quality: PP non-woven fabric

Fabric Type: SS or SMS fabric

Usage: Standard protective coverall

Sterilization: Non-sterile

Weight:

Classification:

Certifications :

Size: Free Size

Colour: Blue

Product features:

Made of polyethylene or polypropylene non-woven 30 g/m2 or 40 g/m2 (for shoe covers with anti-skid layer)

Elastic rubber band at the opening

Protection against contaminating rooms where higher cleanliness

Standards are required

Increased durability

Disposable

Spunbond (SS or SMS) consists of 100% PP non-woven fabric

Antistatic

Dust-proof

Performance indicators: Technical requirements for disposable laminated shoe cover

Application: Hospital, medical

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).

DISPOSABLE POWDER FREE GLOVE





Product Details

Documents: CE, ISO 13485

Brand Name: Dexion Medical

Model Nr: DXNMD-DNG01

Certificate Nr:

Product Name: Disposable Gloves

Use: Disposable gloves. The products are disposable gloves.

Material: % 100 Spunbond PP non-woven shoe cover

Material quality: PE (polyethylene) film-coated non-woven fabric

Sterilization: Non-sterile

Weight:

Classification: Class CP, Class 100, Class10.

Size: XS, S, M, L, XL

Type: Powder-Free

Colour: Blue,White

Product features:

This products are disposable gloves

It is made of synthetic nitrile rubber.

100 % Latex Free.

100 % Synthetic Nitrile

Protein free

No allergy

Performance indicators: Technical requirements for disposable gloves

Application:

The products are non-sterile and are for medical, widely used in fields of medical examination, dentistry, first-aid etc.

Ambidextrous, palm and finger textured, beaded cuff, off line double chlorinated with cleanroom process.

Product Advantages:

More exceeding feature of puncture-resistance, anti-bacteria's penetration, chemical-proof.

Durable & Flexible, Finger-textured, Soft feeling, Comfortable donning.

The tapered cuffs are easier to wear and operate.

Good protective and physical properties, better than latex gloves.





Nontoxic, Harmless and Odorless.

Choiceness Formula, Advanced Technology, Soft Feel, Comfortable, Skid Resistance and Flexible



Syringe List



Syringe List

Product Name	Description	Picture	Spec.
Validate 7days			
Self-destruct syringe (pull and break)	Blister Packing		1~2ML
			2.5~3ML
			5ML
			10 ML
Self-destruct syringe (steel sheet)	Blister Packing		0.05ML
			0.1 ML
			0.5 ML
Validate 7days			
Safety syringe (sliding sleeve)	Blister Packing		1ml
			2~3ml
			5ml
			10ml
Safety syringe (needle cover)	Blister Packing		1ml
			2~3ml
			5ml
			10ml
			20ml

Syringe List

Product Name	Description	Picture	Spec.
Validate 7days			
"disposable syringe"	Blister Packing		1ml
			2~3 ML
			5 ML
			10 ML
			20 ML
			30 ML
			50~60ML
Validate 7days			
disposable needle	Blister Packing		30G
			29G
			28G
			27G
			26G
			25G
			24G
			23G
			22G
			21G
			20G
			18G
16G			

Shield & Glasses



Shield & Glasses



Gallery











Production Place





Yenibosna Merkez Dist.
29 Ekim Street. No:3
IVP Plaza:1 Floor:8
BAHCELİEVLER/ İSTANBUL / TURKEY

info@dexxonmedical.com

www.dexxonmedical.com