



**SMS Surgical Gown,
Disposable Patient Use
SMMS Scrub Suit, Non-Sterile**

FDA **CE** **MADE
IN TURKEY**

Product Details

Product Description : SMS Surgical Gown / Scrub Suits

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSSS14

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21061104

Product Name: Disposable non-sterile short-sleeved patient gown

Type : SMS Surgical Gown / Scrub Suits

Usage: Standard patient gown

Sterilization: Non-sterile

Weight: 35-40 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown / Scrub Suits

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS



SMS Surgical Gown, Disposable Patient Use SMMS Scrub Suit, Non-Sterile

FDA CE MADE
IN TURKEY

Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level : Level 1

Feature:

Disposable and Breathable

Spunbond (SMS) consists of 100% PP non-woven fabric

Economical disposable short sleeves

Antistatic

Each product is packed individually

EN 13795

Style/Pattern: Two-piece set including shirt and pants.

Shirt:

Collar: Around collar or V-collar

Sleeves: Without or short sleeves

Pocket: Available

Pants: With tie on the waist, pockets are available.

Application: Hospital, laboratory etc

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).